



Magnify™ Kit User Guide

FFPE tissues, v1.26.05



Video Demo

Step 1. Follow standard protocol (e.g., <https://www.abcam.com/protocols/ihc-deparaffinization-protocol>) to de-paraffinize the FFPE tissue slide. Use a Kimwipe to remove any residual liquid surrounding the tissue. Choose the smallest gel chamber sticker that accommodates the tissue size and adhere it to the dried surface around the tissue.




Step 2. Take out Solution A and Solution G from the fridge. Pipette 500 μ L of Solution G into an empty 2 mL Eppendorf tube (not provided).


Step 3. Pipette 600 μ L of Solution A into a provided Activation tube and vortex to make the **Activation Buffer**.

Step 4. Pipette 25 μ L of the **Activation Buffer** to the 500 μ L of Solution G to make **Activated Gelling Solution**.


Caution: The Activated Gelling Solution is not stable and it should be used within 15 mins before premature gelling happens.




Step 5. Fill up the gelling chamber with the **Activated Gelling Solution**. Add carefully and try not to perturb the tissue. Place the tissue slide in the fridge or on ice/ ice pad and wait for 15 minutes to allow sufficient infiltration. **Caution: Ensure that the activated gelling solution fills the chamber fully or even slightly above the level of the gel chamber (see video demo). It is crucial to prevent incomplete gelling and the formation of bubbles or gaps within the chamber.**




Step 6. **Gelling:** Place the tissue slide in a Petri dish. Cover the slide with the provided plastic cover. Place the wet Kimwipe and close the lid of the Petri dish to create a humid environment. Incubate the slide at 37°C for 2 hours (preferred), OR 45 °C for 1.5 hours. **Caution: It is critical to ensure the plastic sheet cover not drift as it can introduce bubbles.**



Step 7. **Gel detachment:** Take out the Petri dish from the incubator. For temporary storage before homogenization, store the Petri dish at 4 °C. Remove the plastic sheet cover from the tissue slide. Use a razor blade to cut along the **gel chamber's border** to cut the gel free (see video). Take off the gel chamber sticker and use the blade to clean off any leftover gel outside the chamber. Gently detach the gel with a brush or a gloved finger; it should come off smoothly without dampening it unless it is completely dried.




Step 8. **Homogenization:** Carefully transfer the detached gel into a 15 mL Falcon tube (not provided). Add 3 mL of Buffer H, ensure the gel is fully submerged, and incubate at either 95 °C or autoclave at 120 °C, 12 psi (preferred). Incubation time varies depending on tissue type - consult the tissue lookup table for the specifics. If you see precipitates in Buffer H, please warm them to 45 °C for a few minutes, and they will dissolve quickly. **Caution: It is critical to ensure adequate incubation time and temperature for full homogenization. The times listed in the table are just for reference and It might need further optimization. Samples should appear clear after complete homogenization.**



Tissue Type	Homogenization Time	
	95 °C	120 °C, 12 psi
FFPE Optic Nerve	8 Hours	120 Minutes
FFPE Prostate	8 Hours	90 Minutes
FFPE Brain	100 minutes	Not recommended
All Other FFPE	6 Hours	80 Minutes

Step 9. **Buffer Preparation:** Mix 5 mL of 10X Buffer W with 45 mL of water to create 1x Buffer W solution.

Step 10. **Post-homogenization washing:** Place the specimen in the gel catcher to discard the homogenization buffer. For small gelled tissue, use a six-well plate for easier washing. Wash five times with **1X Buffer W** on an orbital shaker at 70 rpm, 15 minutes each, using 5 mL of fresh 1X Buffer W per wash. After washing, gels can be stored at 4 °C in 1X Buffer W or 1X PBS. **Caution: Sufficient washing is vital to remove surfactants and ensure staining quality. For larger gelled tissue, a 15 mL Falcon tube is recommended to provide more space for effective washing.**




Wash on orbital shaker




Wash 5 times @ RT
15 min/wash
5 mL fresh 1X buffer /wash

Step 11. **Staining:** Prepare 1-2 mL of staining solution with labeling reagents of the user's choice using **Buffer S**. For immunostaining, we recommend diluting the antibody or lectin with buffer S to 1 $\mu\text{g}/\text{mL}$ and diluting lipid stain and pan-protein stain NHS-ester with buffer S to 2 $\mu\text{g}/\text{mL}$ and diluting FISH probes to 10 μM . Incubate the gel with the prepared staining solution in the 15 mL Falcon tube for 1-3 hours at RT. **Caution: The NHS-ester pan-protein stain is known to modify epitopes; therefore, it is advised to apply this stain to the specimen only after all other staining procedures have been completed.**


Step 12. **Post-staining washing:** Wash the specimen with 1X Buffer W or 1X PBS again three times, like in step 10. Optionally, User can stain the specimen with DAPI solution.



Step 13. **Expansion:** Transfer the gel to the one-well glass bottom plate, ensuring the orientation of the cut corners as illustrated to keep the objective lens close to the tissue during imaging. If the gel is incorrectly oriented, gently correct it by adding 1 mL of 1X PBS and using the brush to flip it carefully. It's crucial to handle the gel delicately to avoid damage. Fill up the plate with 1X PBS for 4X expansion on the orbital shaker for 15 mins. Use droppers to carefully remove all liquid before imaging.



Step 14. For 11X expansion, wash the gel inside the one-well plate with fresh de-ionized water using droppers on the orbital shaker at least three times for 15 minutes each. After the gel has fully expanded, use droppers to carefully remove all liquid before imaging.



Step 15. To minimize drifting during imaging, user may optionally coat the sample with a warm 1% agarose solution in ddH₂O, ensuring to envelop both the sample's surroundings and its top. Wait for the agarose to solidify fully before starting the imaging process.

Disclaimer

Restrictions and Liabilities

This technical data sheet is provided “as is” and Magnify™ Biosciences assumes no responsibility for any typographical, technical, or other inaccuracies. The document is subject to change, without notice, in future editions.

To the maximum extent permitted by applicable law, Magnify™ Biosciences disclaims all warranties, either express or implied, with regard to this technical data sheet and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Magnify™ Biosciences shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

This document might contain references to third party sources of information, hardware or software, products, or services and/or third party web sites (collectively the “Third-Party Information”). Magnify™ Biosciences does not control, and is not responsible for, any Third-Party Information. The inclusion of Third-Party Information in this document does not imply endorsement by Magnify™ Biosciences of the Third-Party Information or the third party in any way.

Magnify™ Biosciences is not responsible nor will be liable in any way for your use of any software or equipment that is not supplied by Magnify™ Biosciences in connection with your use of Magnify™ Biosciences products.

Magnify™ Biosciences does not in any way guarantee or represent that you will obtain satisfactory results from using Magnify™ Biosciences products as described herein. The only warranties provided to you are included in the Limited Warranty enclosed with this document. You assume all risk in connection with your use of Magnify™ Biosciences products.

Note to Purchaser: Limited Product Warranty

Any product that does not meet the performance standards stated in the product specification sheet will be replaced at no charge. This warranty limits our liability to the replacement of the product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Magnify™ Biosciences. Magnify™ Biosciences shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use any product.

Note to Purchaser: Limited License

Magnify™ Expansion Kits are developed, designed, and sold exclusively for research purposes. Neither the product, nor any individual component, has been tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to the SDS, which is available on request.

Certain applications of this product are covered by patents issued to parties other than Magnify™ Biosciences and applicable in certain countries. Purchase of this product does not include a license to perform any such applications. Users of this product may therefore be required to obtain a patent license depending upon the particular application and country in which the product is used.